



CAVALLO s.r.l.
via Mulino 5/B
20090 BUCCINASCO (Milano) ITALY
T.+39.0245101515 Fax +39.0245100234
e-mail: info@cavallocom.it www.cavallocom.it



REFRIGERATED INCUBATOR model 1400FTX 4T with n.4 Independent Temperature Chambers



Technical features

Temperature range: **+5°C to +50°C**

Outer case in stainless steel as per AISI304.

High density polyurethane insulation.

n.4 Internal independent Chambers in stainless steel as per AISI 304

n.4 External front Doors in stainless steel as per AISI 304

CFC-free refrigerator unit for ozone and environment protection (Art.549/93)

Fan-forced air circulation

Heating provided by stainless steel sheathed resistors

Inner dimensions of each Chamber: cm.57x69x69 (wxdxh)

net capacity of each chamber: 271 lt.

N°2 stainless steel Shelves (total shelves:8)

Control Panel housing:

n.4 Digital temperature microprocessors P.I.D. auto-tuning control

n.4 over-temperature safety thermostat with alarm light

n.4 Mains Switch

Safety fuses

230V/1 50 Hz Electric system acc. to CEI 62-25 - CE Marking

Model 1400FTX 4T

External dimensions: cm 140x85x206 (wxdxh)

n.4 casters

Optional:

Ø 35mm hole/port

Made in Italy