



CAVALLO s.r.l.
via Mulino 5/B
20090 BUCCINASCO (Milano) ITALY
T.++39.0245101515 Fax +39.0245100234
e-mail: info@cavallocom.it www.cavallocom.it



STERILIZING DRYING OVEN



Technical features:

Temperature range: from Ambient +5°C to +250°C
Outer case in stainless epoxy painted steel plate with acid and rust-proof ground
Internal chamber in stainless steel as per AISI 304
External front door in enamelled steel and interior in stainless steel with silicone rubber gasket
Thermo insulation through thick rock wool panels
Heating provided by n.2 stainless steel sheathed resistors easy to replace
Fan powered circulation, through n.2 specific fans suitable for high temperature, ensure excellent thermal uniformity
Exhausting duct with movable cover on top
Control Panel mounted on the right side, housing :
-Digital thermoregulator PID autotuning
-Pt100 Sensor
-Over-temperature safety thermostat with alarm light
-Mains Switch
-Alarm light
-Safety fuses
230/1/50 Hz Electric system according to CEI 62-25 – CE Marking

Mod. C4

int. dimensions: cm. 64x65x100H = lt.416
ext.dimensions: cm.102x79x123H
n.3 st.steel Shelves

Options:

Timer
Key Lock
Castors kit

Mod. C5

int. dimensions: cm. 64x65x140H = lt.582
ext.dimensions: cm.102x79x164H
n.4 st.steel Shelves

Drying oven version
IQ-OQ Protocols - Temperature Protocol -
Temperature Recorder
Customized version upon request

MADE IN ITALY