



CAVALLO s.r.l.
 via Mulino 5/B
 20090 BUCCINASCO (MI)
 T.+39.02.45101515 F.+39.0245100234
 e-mail:info@cavallocom.it www.cavallocom.it



GLASSWARE DRYER



Technical features:

- Temperature range: Ambient +5°C to 150°C
- Outer case in stove enamelled steel plate
- Internal chamber in stainless steel case 18/8 AISI 304
- Front door with 2 shutters, with interior in stainless steel and exterior in stove enamelled steel plate
- Fan-forced air circulation via centrifugal fan
- Outlet for the saturated air through a hole in the body upper part
- Heating via In Cr cable
- Thermo-insulation through thick rock wool panels
- Internal shelves in stainless steel
- Side Control Panel housing:
 - Digital thermostat with microprocessor allowing P.I.D. auto-tuning control
 - Safety thermostat with visual alarm
 - Digital timer for cycle setting
 - Timer switch
 - Line switch
 - Safety fuses

V.380/3 +N.+Ground (or V.230/1)
 CEI 62-25 Electric device – CE Marking

Model	internal dimensions (wxdxh) cm	total Volume ltrs	Shelves (tot.capacity)	external dimensions (wxdxh) cm	Watt
3671 A	90x60x150	810	4 (12)	130x 80x219	4500
3671 B	110x70x150	1155	4 (12)	150x 90x219	6000
3671 C	130x80x150	1560	4 (12)	170x100x219	7000

MADE IN ITALY

Data, description and illustration are only indicative, products may have different characteristics due to manufacturing requirements. Cavallo srl reserves the right to any modifications to its products without notice.