



CAVALLO s.r.l.
via Mulino 5/B
20090 BUCCINASCO (Milano) ITALY
T.++39.0245101515 Fax +39.0245100234
e-mail: info@cavallocom.it www.cavallocom.it



Refrigerated Incubator model CTA



Technical features:

Temperature range: +5°C/+50°C
Outer case in stainless epoxy painted steel plate with acid and rust-proof ground
Internal chamber in stainless steel as per AISI 304
External front door in enamelled steel and interior in stainless steel with silicone rubber gasket
Internal Glass door
Thermo insulation through thick rock wool panels
Heating through stainless steel sheathed Resistor easy to replace
CFC-free refrigerator unit for ozone and environment protection (Art.549/93)
Overpowered compressor
Fan powered circulation, through specific fan to ensure excellent thermal uniformity
Control Panel mounted on the right side, housing :
-Ascon PID autotuning microprocessor Regulator
-Under and over-temperature safety thermostats with alarm light -Mains Switch
-Alarm light
-Safety fuses
n.2 Stainless steel Shelves

230/1/50 Hz Electric system according to CEI 62-25 CE Marking

Mod. C2TA Mod C3TA

int. dimensions: cm. 46x46x60 H = lt.125 int. dimensions: cm. 61x63x60 H = lt.230
ext.dimensions: cm. 85x62x76 H ext.dimensions: cm.102x80x76 H

Options:

IQ-OQ Protocols
Temperature Protocol
Temperature Recorder
Remote Alarm system for Temperature through NO/NC 'clean contact'
Steel Stand with castors

Made in Italy