



CAVALLO s.r.l.
via Mulino 5/B
20090 BUCCINASCO (MI)
T.+39.02.45101515 F.+39.0245100234
e-mail:info@cavallocom.it www.cavallocom.it



Air INCUBATOR



Technical features:

Temperature range: from Ambient +5°C to +60°C

Outer case in stainless epoxy painted steel plate with acid and rust-proof ground

Internal chamber in stainless steel 18/8 as per AISI 304

External front door in enamelled steel and interior in stainless steel with silicone rubber gasket
inner transparent Door

Thermo insulation through thick rock wool panels

Heating provided by n.2 stainless steel sheathed resistors easy to replace

Exhausting duct with movable cover on top

Control Panel mounted on the right side, housing :

- Digital Microprocessor Regulator PID autotuning
- Over-temperature safety thermostat with alarm light
- Mains Switch
- Alarm light
- Safety fuses

230/1/50 Hz Electric system according to CEI 62-25 – CE Marking

Mod. C4TE

int. dimensions: cm. 64x63x90H = lt.360

ext.dimensions: cm.102x80x123H

n.3 st.steel Shelves

Mod. C5TE

int. dimensions: cm. 64x63x130H = lt.522

ext.dimensions: cm.102x80x164H

n.4 st.steel Shelves

Options:

Key Lock

Castors kit

ventilated/fan circulated version

IQ-OQ Protocols - Temperature Protocol - Temperature Recorder

Customized version upon request

MADE IN ITALY

Data, description and illustration are only indicative, products may have different characteristics due to manufacturing requirements.
Cavallo srl reserves the right to any modifications to its products without notice.