



**CAVALLO s.r.l.**  
via Mulino 5/B  
20090 BUCCINASCO (Milano) ITALY  
T.++39.0245101515 Fax +39.0245100234  
e-mail: [info@cavallocom.it](mailto:info@cavallocom.it) [www.cavallocom.it](http://www.cavallocom.it)



## Shaking Waterbath



### Technical features:

24/24 HOUR FUNCTION GUARANTEED

Outer case in stainless epoxy painted steel plate with acid rust-proof ground

Bath made in stainless steel 18/8 as per AISI 304

Glass-wool thermo-insulated walls and bottom

Inner movable trolley mounted on noiseless auto-lubricant wheels

Motion deriving from a motor combined with gear reduction motor

Shaking speed between 30/180 strokes per minute

Keyboard for speed setting

LCD display for direct reading of the number of stroke per minute

Heating via sheathed immersion resistors (1400 Watt )

Rack holder in stainless steel

Front control panel housing:

digital thermostat with microprocessor allowing P.I.D. auto-tuning control

safety thermostat with indicator lamp

shaking speed electronic controller and read-out

heating switch

shaking switch

Ext. dimensions (WxDxH) 870x380x350 mm = 36 kgs

Bath dimensions (WxDxH) 550x320x230 mm = 40 ltrs

Rack dimensions 400x270

CEI 62-25 Electric device – CE Marking - V.230 singlephase 50 Hz - 1500 Watt

**model DUBNOFF** (crosswise shaking motion)

**model ROTABATH** (orbital shaking motion)

3161 - Rack for 40 pcs. E-flasks 25 ml

3162 - Rack for 35 pcs. E-flasks 50 ml

3163 - Rack for 24 pcs. E-flasks 100 ml

3164 - Rack for 12 pcs. E-flasks 250 ml

3165 - Rack for 6 pcs. E-flasks 500 ml

3166 - Rack for test-tubes  $\varnothing$  on request

3167A - Lid

3167B - Cooling Coil

CV1 - Steel plate Stand

**MADE IN ITALY**