



CAVALLO s.r.l.
via Mulino 5/B
20090 BUCCINASCO (Milano) ITALY
T.+39.0245101515 Fax +39.0245100234
e-mail: cavallo@cavallocom.it www.cavallocom.it



REFRIGERATOR/FREEZER with two separated cabinets



Refrigerator/Freezer

with two separated cabinets
air ventilated

model CA700TN/BT NEW

2 solid Doors

2 traditional motors (multiple temperature)

Temperature range: -2°C /+8°C

Temperature range: -15°C /-22°C

Structure made entirely of stainless steel as per AISI 304 with Scotch-brite exterior finish.

Rounded inner corners for easy cleaning operations, bottom bowl to contain spills.

Stainless steel Door. easy cleaning magnetic removable gasket. Self-closing system for openings <90°.

Ergonomic Door Handle. Positive seal self closing Door. Lock and key.

The insulation thickness is 60 mm polyurethane injected CFC and HCFC.

Internal service lighting: Led light activated when door open.

Traditional Refrigeration unit: continuous forced ventilation, located at the top of the structure.

The ventilated refrigeration system allows a good distribution of set temperature inside the chamber.

Defrosting: automatic hot gas

Double temperature probe.

Control panel: Touch screen front panel. Keypad lock function. Door open alarm. Possibility to select temperatures in °C or °F.

Fast overcooling function.

Compressor delivered cooling: Watt 250/287

Absorption: Watt 304/362

Refrig. gas: R404a

n.4 plastic-coated Shelves 650x530mm. with adjustable rails.

n.4 adjustable feet

External dimensions mm.: 700x850x2060 LxPxH = Capacity: 700 ltrs

net weight Kgs: 140

Electricity (V/ph/Hz): 230/1/50HZ – CE marking

Made in Italy