



**CAVALLO s.r.l.**  
via Mulino 5/B  
20090 BUCCINASCO (Milano) ITALY  
T.++39.0245101515 Fax +39.0245100234  
e-mail: [info@cavallocom.it](mailto:info@cavallocom.it) [www.cavallocom.it](http://www.cavallocom.it)



## STERILIZING DRYING OVEN



### Technical features:

Temperature range from Ambient +5°C to +250°C  
Outer case in stainless epoxy painted steel plate with acid and rust-proof ground  
stainless steel Internal chamber 18/8 AISI 304  
External front door in enamelled steel and interior in stainless steel with silicone rubber gasket  
Thermo insulation through thick rock wool panels  
Heating through stainless steel sheathed Resistor easy to replace  
Heating provided by stainless steel sheathed resistors  
Fan powered circulation, through specific fan suitable for high temperature, ensures excellent thermal uniformity  
Exhausting duct with movable cover on top  
Control Panel mounted on the right side, housing :  
-Digital thermoregulator PID autotuning  
-Pt100 Sensor  
-Over-temperature safety thermostat with alarm light  
-Mains Switch

-Alarm light  
-Safety fuses  
230/1/50 Hz Electric system according to CEI 62-25 – CE Marking

### Mod. C2

int. dimensions: cm. 47x48x60H =  
lt.135  
ext.dimensions: cm. 85x62x75H  
n.3 grid Shelves

### Mod. C3

int. dimensions: cm. 64x65x60H =  
lt.249  
ext.dimensions: cm.102x78x75H  
n.2 st.steel Shelves

### Options:

Timer  
Key Lock  
Steel Stand with castors  
IQ-OQ Protocols - Temperature Protocol -  
Temperature Recorder  
Customized version upon request

**MADE IN ITALY**