



CAVALLO s.r.l.
via Mulino 5/B
20090 BUCCINASCO (Milano) ITALY
T.++39.0245101515 Fax +39.0245100234
e-mail: info@cavallocom.it www.cavallocom.it



STERILISING OVEN



Double door oven, specifically designed for pharmaceutical labs, suitable to be accommodated in places between sterile and pre-sterile rooms.

Technical features:

- Temperature range above ambient + 5°C up to 250°C
- Outer case and internal chamber in stainless steel according to AISI 304
- Stainless steel grid shelves, ladder-rack type, to improve air flow and provide an even more temperature distribution within the oven
- Glass fibre insulation
- Double door featuring chromate brass hinges and closing onto silicone rubber gaskets
- Manual or automatic door lock with visual indicator lamp
- Double chamber for forced circulation of hot air via electric fan suitable for high temperature
- Heating via battery of flyer resistors in stainless steel

Side Control Panel housing:

- Digital thermostat with microprocessor allowing P.I.D. auto-tuning control
- Alarm light for cycle ending
- Over-temperature safety thermostat
- Multi-scale timer
- Line switch
- Heater switch
- Safety fuses
- Acoustic and visual alarm light signalling cycle ending

380 V + N + Ground

Different dimensions are available upon request

Electric supply according to CEI 62-25